ET KASUTUSJUHEND

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GB INSTRUCTIONS FOR USE

General notes

Alcotest can be used only once. Alcotest is designed to measure breath alcohol concentration. Alcotest detects blood alcohol concentration of 0,2%*, i.e. breath alcohol concentration of 0,1 mg/l. The manufacturer shall not be liable for the result or reliability of the test. The alcohol level increases up to approximately an hour from alcohol consumption. Even less than 0,2%* of alcohol in your blood may impair your ability to drive, causing an immediate threat to your life and health and a risk to road users.

The alcohol concentration reading may not be used as evidence, particularly as evidence confirming being/not being in a state of intoxication/under the influence of alcohol in any proceedings. The Alcotest user shall be solely responsible for the interpretation of the test results. The test results shall not be used to question the results of tests performed by officers or at their request. The manufacturer shall not be liable for using the product not in line with these Instructions for Use.

Components

Alcotest consists of a glass tube, measuring substance – reagent, adsorption granules, filters, aluminium foil and two protective caps. Individual Alcotest components may pose a threat to life or health. When used incorrectly, the glass tube may break and cause cuts. Other elements may be swallowed. Throw the used Alcotest into a municipal waste bin. The used Alcotest poses no threat to the environment according to the European Union Directive. Do not leave the used Alcotest in a place, where it may be taken away by an unauthorised person.

Steps:

1.Push the protective caps on both ends of Alcotest by pressing strongly./Pierce mouthpiece both sides of the breathalyzer.

2.Depending on the Alcotest model, insert the tube into the balloon. If there is no balloon, move on to the next step.

3.Take a deep breath and blow the air into the glass tube (twice, for about 10 seconds). If a disposable breathalyser includes a balloon, blow the air into the glass tube until the balloon is fully inflated.

4.Wait at least two minutes. If the colour in the central part of the tube changes, it means that you are under the influence of alcohol. After 5 minutes, the result may be unreliable.

5.Compare the colour with the colour scale on the packaging. The most accurate comparison results are obtained in daylight. The Alcotest reading similar to the one marked 0,2%* on the packaging indicates that the blood alcohol concentration is 0,2%* or more.

Precautions

-Keep Alcotest out of reach of children.

-It is forbidden to remove the measuring substance from Alcotest.

-It is recommended not to pour a windshield washer fluid containing alcohol before using Alcotest, because the liquid vapours may affect the reading.

-Avoid direct contact of the measuring substance with eyes and skin.

-Do not ingest the measuring substance.

-Do not throw Alcotest into fire.

-In the case of contact of the measuring substance with eyes, rinse the eyes with plenty of water while keeping the eyelids wide open. Immediately seek medical advice. It is also recommended to consult an ophthalmologist.

-In the case of contact of the measuring substance with skin, rinse the affected area with plenty of water. In the event of any skin reactions (redness, rash, allergic reactions or wounds), immediately seek medical advice. It is also recommended to consult a dermatologist.

-If measuring substance is ingested, drink plenty of water. Do not induce vomiting. Immediately seek medical advice.

-If the measuring substance is inhaled, go to fresh air. Seek medical advice if needed.

Contains sulfuric acid and phosphoric acid.

-Protect from sunlight and high temperatures.

-Store in a dry and safe place out of reach of children.